

genomica.com


# GENOMICA

MULTIPLEX · MICROARRAY MDX



CE

VERSION 3 ■ APRIL 2014



INNOVATION  
& RELIABILITY

## CLART<sup>®</sup> by GENOMICA

GENOMICA's Clinical Array Technology, CLART<sup>®</sup>, is a diagnostics platform based on a Multiplex-PCR amplification step and the subsequent analysis of the amplified material on a low density microarray (LD-Microarray).

The combination of both allows the simultaneous detection of multiple targets in a single assay, providing detailed information and thus facilitating the decision-making process for the clinicians.

The analysis of the sample is simple, straightforward and robust making it suitable for every molecular diagnostic laboratory.



## CAR<sup>®</sup>

Clinical Array Reader is a colorimetric array reader unit, running our proprietary software, SAICLART<sup>®</sup>, for the analysis and interpretation of the array images. SAICLART<sup>®</sup> has been designed and validated for interpreting the arrays images turning them automatically into clinically relevant data.

The reader displays an interface called CLEIS (CLART<sup>®</sup> End-user Interface Software), based on the extensive experience and the customers feedback provided, thus obtaining a very intuitive user-friendly format.

### Technical characteristics:

- > Integrated PC.
- > Touch screen.
- > Easy data management:
  - LIS bidirectional connection.
  - HTML and bmp formats.
  - Printable, exportable and storable reports.




## SAICLART<sup>®</sup>

SAICLART<sup>®</sup> is GENOMICA's proprietary image processing software for the analysis and interpretation of microarrays images. Its advanced recognition algorithm is able to discern spots from artifacts like speckles, fibers or bubbles, greatly reducing false positive and false negative results.

The image analysis is completely automated, avoiding subjectivity as result of user interpretation, thus providing fast, accurate and reproducible results, improving the overall performance of the assay.

## autoclart<sup>®</sup>

Autoclart is the latest GENOMICA's automation system. It's designed for microplate processing and provides reagents addition, liquid handling, heating, cooling and microplate shaking.



Its technical features makes it suitable for almost all molecular diagnostics laboratories throughput:

- > Compact size: 60 x 60 x 65 cm.
- > Plug and play system.
- > Operator interface is via an integral control panel containing a knob and a LCD display screen.
- > Microplate holder offers: heating and cooling under a peltier controlled system, giving precise temperature control during assay, but also shaking wells during incubations.

### Autoclart main advantages:

- > Notably reduces the hands-on time required for performing any CLART<sup>®</sup> assay.
- > Minimizes the potential intra-assay variability, but also inter-laboratory variability.
- > All consumables are housed inside the instrument, giving a very small footprint in the lab.
- > Flexibility in number of samples per run, from 4 to 96 at once.



# CLART® CMA

## KRAS·BRAF·PI3K & NRAS·iKRAS

Specific detection of somatic mutations in oncogenes determining response therapy

### MUTATIONS DETECTED



### MAIN ADVANTAGES OF MULTIPLE DETECTION OF SOMATIC MUTATIONS

- Detects 20 of the most prevalent mutations on KRAS, BRAF, PI3K, NRAS and iKRAS.
- Reduces the amount of sample required
- All mutations included in one kit can be detected in a single array.
- Avoids unnecessary toxicity caused by improper selected antitumor therapy, as well as its associated costs.

### FEATURES


- Both kits have been validated for automatic and manual DNA extraction from FFPE samples and cell lines.
- Mutational status can be detected for single or multiple genes just combining references.
- High sensitivity and specificity.
- Three internal quality controls included per sample:
  - **Genomic DNA control:** validates the extraction performance.
  - **Amplification control:** avoids false negative results.
  - **Biotin markers:** check the proper performance of the visualization reagents provided with the kit.
- Each mutation is detected in triplicate avoiding unspecific bindings.
- Results are obtained within a working day.

### Data Management

- Automatic reading and interpretation of results (CAR®).
- User-friendly report format (html, bmp)
- Printable, exportable and storable files.

### REPORTING RESULTS

> Report and image obtained by CAR® reader.



**A1**

KRAS Analysis - Not tested			
Controls	Result	Controls	Result
Extraction	Not tested	Amplification	Not tested
Mutations			
Mutation	Result	Mutation	Result
KRAS 12 Ser	Not tested	KRAS 12 Arg	Not tested
KRAS 12 Asp	Not tested	KRAS 12 Cys	Not tested
KRAS 13 Asp	Not tested	KRAS 61 Leu	Not tested
KRAS 61 His (CAT)	Not tested	KRAS 12 Val	Not tested
KRAS 12 Ala	Not tested		

**BRAF Analysis - Not tested**

### ORDERING REFERENCES

#### CLART® CMA KRAS-BRAF-PI3K

CLART® CMA KRAS  
Amplification 24 tests: CS-0412-24

CLART® CMA BRAF  
Amplification 24 tests: CS-0512-24

CLART® CMA PI3K  
Amplification 24 tests: CS-0612-24

CLART® CMA KBP Array  
Visualization 24 tests: CS-0712-24

#### CLART® CMA NRAS-iKRAS

CLART® CMA NRAS-iKRAS  
Amplification 24 tests: CS-0114-24

CLART® CMA NiK Array  
Visualization 24 tests: CS-0214-24


- > CLART® CMA KRAS-BRAF-PI3K: Essential features of CLART® CMA KRAS-BRAF-PI3K are protected by European Patent Application N° EP11382397.5 and EP12382370.0
- > CLART® CMA KRAS-BRAF-PI3K and CLART® CMA NRAS-iKRAS kits accomplish with all the normative described in EU 98/79/EC for IVD.
- > CLART® CMA KRAS-BRAF-PI3K is included in the recommended methods for testing BRAF V600E mutation by the Danish Pathology Association\*. <http://danskspatologi.dk/?pid=191>



# CLART® HPV2

## Genotyping Human Papillomavirus

### GENOTYPES DETECTED



> High risk & low risk distribution detected by GENOMICA according to Dunne, E. JAMA 2007;297(8):813-81.

### IMPORTANCE OF HPV GENOTYPING

- Genotyping allows simultaneous detection of single infections or co-infections.
- Provides information about HPV prevalences, specially among already vaccinated cohorts.
- Enables the early detection and patient follow up, essential for cancer prevention.
- Allows studies about HPV types and their distribution in rectal, pharyngeal and cervical cancer.

### FEATURES

- Is the first genotyping assay used in HPV screening programs worldwide thanks to GENOMICA's automated system (autoclart®).
- Detects and genotypes 35 different HPV types, including HR and LR in one assay.
- Kits are validated for automatic and manual extraction of LBC, Swabs and FFPE tissues.
- High sensitivity and specificity. Clinical validation performed.



- Three quality controls included per sample:
  - **Genomic DNA control:** validates the extraction performance.
  - **Amplification control:** avoids false negative results.
  - **Biotin markers:** check the proper performance of the visualization reagents provided with the kit.
- Each HPV genotype is detected in triplicate avoiding unspecific bindings.
- Results are obtained within a working day.

### Data Management

- Automatic reading and interpretation of results (CAR®).
- User-friendly report format (html, bmp).
- Samples are processed individually and three complementary reports are generated.
- Printable, exportable and storable reports.

### REPORTING RESULTS

> Report and image obtained by CAR® reader.



C3		
Result view <b>CLART® HPV2</b>		AT code Rev. 06/13
Sample reference:	3	
Assay ID:	00000050516 (C1)	
Analysis type:	tmf end point detection	
Date and time:	Fri, Jan 14 15:25:36 2011	
Virus	Result	Controls
Type 6	Negative	Passed
Type 11	Negative	Passed
Type 16	Negative	Passed
Type 18	Positive	Passed

### ORDERING REFERENCES

CLART® HPV2 Extraction

48 tests: AT-1105-48

CLART® HPV2 Amplification

48 tests: AT-1106-48-MT

CLART® HPV2 Visualization

48 tests: CS-0208-48

> CLART® HPV2: Essential features of CLART® HPV2 are protected by Patent Families of International PCT Patent applications WO2007017699 and WO2011116797.

> CLART® HPV2 kit, accomplish with all the normatives described in EU 98/79/EC for IVD.



# CLART® STIs - A&B

## Detecting microorganisms causing urogenital tract infections

### MICROORGANISMS DETECTED

<b>BACTERIA</b> <ul style="list-style-type: none"> <li>Chlamydia trachomatis*</li> <li>Neisseria gonorrhoeae*</li> <li>Mycoplasma genitalium*</li> <li>Mycoplasma hominis</li> <li>Ureaplasma parvum</li> <li>Ureaplasma urealyticum</li> <li>Treponema pallidum</li> <li>Haemophilus ducreyii</li> </ul>	<b>VIRUSES</b> <ul style="list-style-type: none"> <li>HSV1</li> <li>HSV2</li> </ul>	<b>PARASITE</b> <ul style="list-style-type: none"> <li>Trichomonas vaginalis*</li> </ul>	<b>FUNGI</b> <ul style="list-style-type: none"> <li>Candida albicans</li> <li>Candida glabrata</li> <li>Candida krusei</li> <li>Candida dubliniensis</li> <li>Candida guilliermondii</li> <li>Candida parapsilosis</li> <li>Candida tropicalis</li> </ul>
---	---	--	---

\* Microorganism analyzed with CLART®STI A, remaining microorganism are analyzed with CLART®STI B

### MAIN ADVANTAGES OF MOLECULAR DETECTION OF MICROORGANISMS CAUSING STI

Molecular diagnostic techniques provide greater sensitivity and standardize the multiple methodologies used so far for the detection of those microorganism causing urogenital tract infections.

Moreover, molecular diagnostics techniques reduce the drawbacks and limitations of conventional detection methods as:

- Low sensitivity shown by cultures.
- Antibody titers variations due to antiviral treatment.
- Difficulties in microscopy detection.

### FEATURES

- Both kits have been validated for automatic DNA extraction from urine samples and swabs (vaginal, cervical, endocervical, urethral and rectal).
- High sensitivity and specificity.


- No previous culture required.
- Three quality controls included per sample:
  - Genomic DNA control:** validates the extraction performance.
  - Amplification control:** avoids false negative results.
  - Biotin markers:** check the proper performance of the visualization reagents provided with the kit.
- Each target is detected in triplicate avoiding unspecific bindings.
- Results are obtained within a working day.
- Considerably reduction of turnaround time allowing the most effective therapy adjustment in the short term.

### Data Management

- Automatic reading and interpretation of results (CAR®).
- User-friendly report format (html, bmp).
- Samples are analyzed individually and three complementary reports are generated.
- Each report generated can be stored in the reader, exported or printed.

### REPORTING RESULTS

> Report and image obtained by CAR® reader.



B4		AT code: 00516	Rev: 06/13
Result view: CLART® STIs A			
Sample reference:	3		
Array ID:	00000050516 (C1)		
Analysis type:	trib end point detection		
Date and time:	Fri, Jan 14 15:25:36 2011		
Bacteria	Result	Controls	
Chlamydia trachomatis	Negative	Passed	
Neisseria gonorhoea	Positive	Passed	
Treponema pallidum	Negative	Passed	
Mycoplasma genitalium	Negative	Passed	

### ORDERING REFERENCES

CLART® STIs Amplification A

48 tests: CS-1112-48

CLART® STIs Amplification B

48 tests: CS-0213-48

CLART® STIs A / CLART® STIs B Visualization

48 tests: CS-1212-48

> CLART® STIs A&B for detection of microorganisms that cause infections of the urogenital tract fulfills the European Directive 98/79/EC for MD products.

> CLART® STIs A: The analysis of Chlamydia trachomatis has been assented by the NB0318.



# CLART® PneumoVir

Detection of respiratory viruses

## VIRUSES DETECTED

Adenovirus	Bocavirus
Metapneumovirus A	Coronavirus 229E
Metapneumovirus B	Enterovirus
Parainfluenza 1	Influenza A, subtyping: Seasonal H1N1 Seasonal H3N2 New H1N1
Parainfluenza 2	
Parainfluenza 3	
Parainfluenza 4 subtyping: Parainfluenza 4 A Parainfluenza 4 B	
Rhinovirus	
RSV A	Influenza B
RSV B	Influenza C

## MAIN ADVANTAGES OF RESPIRATORY VIRUS DETECTION

- Co-infections of several types and subtypes can be detected in the same assay.
- Avoids unnecessary treatments and long hospitalizations.
- Allows virus prevalence studies.

## FEATURES

- Types an Influenza complete panel including: seasonal Influenza A H1N1 and H3N2, generic Influenza A and New Influenza A H1N1.
- The kit is validated for both, automatic and manual DNA/RNA extraction, of BAL, Nasopharyngeal washes and Nasopharyngeal swabs.
- High sensitivity and specificity.


- Two quality controls included per sample:
  - **Amplification control:** avoids false negative results.
  - **Biotin markers:** check the proper performance of the visualization reagents provided with the kit.
- Each virus type is detected in triplicate avoiding unspecific bindings.
- Results are obtained within a working day.

### Data Management

- Automatic reading and interpretation of results (CAR®).
- User-friendly report format (html, bmp).
- Samples are processed individually and three complementary reports are generated.
- Each report generated can be stored in the reader, exported or printed.

## REPORTING RESULTS

> Report and image obtained by CAR® reader.



A2		
Result view		AT code: 0016
CLART® PneumoVir		Rev: 10a.3
Sample reference:	3	
Array ID:	0000005016 (C1)	
Analysis type:	trib end point detection	
Date and time:	Fri, Jan 14 15:25:36 2011	
Virus		
Virus	Result	Controls
Adenovirus	Negative	Passed
Bocavirus (hBoV)	Negative	Passed
Coronavirus Type 229	Negative	Passed
Enterovirus	Negative	Passed

## ORDERING REFERENCES

CLART® PneumoVir Extraction

48 tests: AT-0507-48

CLART® PneumoVir Amplification

48 tests: AT-0607-48-MT

CLART® PneumoVir Visualization

48 tests: CS-0408-48


> CLART® PneumoVir: Essential features of CLART® PneumoVir are protected by Patent Families of International PCT Patent application WO2009144497.

> CLART® PneumoVir kit, accomplish with all the normatives described in EU 98/79/CE for MD.



# PneumoCLART bacteria®

## Detecting bacteria causing respiratory tract infections



### MAIN ADVANTAGES OF MULTIPLE MOLECULAR DETECTION OF RESPIRATORY BACTERIA

- Co-infections with virus and bacteria represents more than 60% of all community acquired pneumonia.
- Multiplex MDx techniques allows the detection of such co-infections, providing greater sensitivity and standardizing the multiple methodologies used so far for the detection of those microorganism causing respiratory tract infections.
- Moreover, Multiplex MDx techniques:
  - Reduce the drawbacks and limitations of conventional detection methods.
  - Allows pathogen-directed treatment.
  - Anticipate illness potential complications.

### FEATURES



- PneumoCLART bacteria® has been validated for automatic DNA extraction from sputum, nasopharyngeal lavages/ exudates/aspirates, BAL and bronchial suction.
- Antibiotic resistance detection (Mec A genes).
- High sensitivity and specificity.
- Three quality controls included per sample:
  - Genomic DNA control:** validates the extraction performance.
  - Amplification control:** avoids false negative results.
  - Biotin markers:** check the proper performance of the visualization reagents provided.
- Each target is detected in triplicate avoiding unspecific bindings.
- Results are obtained within a working day.
- Considerably reduction of turnaround time.

### Data Management

- Automatic reading and interpretation of results (CAR®).
- User-friendly report format (html, bmp).

### REPORTING RESULTS

> Report and image obtained by CAR® reader.



A2		AT code: 6059
Result view	PneumoCLART bacteria	
Sample reference:	3	Rev. Min.3
Array ID:	00000050516 (C1)	
Analysis type:	trib end point detection	
Date and time:	Fri, Jan 14 15:25:36 2011	
Bacteria	Result	Controls
Staphylococcus aureus	Negative	Passed
Haemophilus influenzae	Negative	Passed
Streptococcus pneumoniae	Negative	Passed
Moraxella catharralis	Negative	Passed

### ORDERING REFERENCES

PneumoCLART bacteria® Amplification

48 tests: CS-1013-48

PneumoCLART bacteria® Visualization

48 tests: CS-1213-48

> PneumoCLART bacteria® accomplish with all the normatives described in EU98/79/EC for IVD.





# CLART® ENTERPEX

## Detection of Human Herpesvirus and Enterovirus

### VIRUSES DETECTED

Herpes Simplex Virus 1 (HSV1)
Herpes Simplex Virus 2 (HSV2)
Varicella Zoster Virus (VZV)
Epstein-Barr Virus (EBV)
Cytomegalovirus (CMV)
Human Herpes Virus 6 (HHV6)
Human Herpes Virus 7 (HHV7)
Human Herpes Virus 8 (HHV8)
Enterovirus (Coxsackievirus, Poliovirus and Enterovirus)

### MAIN ADVANTAGES OF MULTIPLEX DETECTION

- Simultaneous detection of those virus causing a wide spectrum of diseases and infections.
- Single infections and co-infections can be detected.
- Best way to detect specific Enterovirus causing diseases, differentiating them from those sharing clinical and epidemiological characteristics.
- Reduced sample volume consuming obtaining the most accurate information.

### FEATURES


- The kit is validated for automatic and manual DNA/RNA extraction of CSF, Serum, Plasma, Swabs and FFPE samples.
- High sensitivity and specificity.
- Two quality controls included per sample:
  - **Amplification control:** avoids false negative results.
  - **Biotin markers:** check the proper performance of the visualization reagents provided with the kit.
- Each virus is detected in quadruplicate avoiding unspecific bindings.
- Results are obtained within a working day.

### Data Management

- Automatic reading and interpretation of results (CAR®).
- User friendly report format (html, bmp).
- Samples are analyzed individually and three complementary reports are generated.
- Each report generated can be stored in the reader, exported or printed.

### REPORTING RESULTS

> Report and image obtained by CAR® reader.



A5		
Result view	A7 min: 00:00	
CLART® ENTERPEX		
Sample reference:	3	
Array ID:	0000005016 (C1)	
Analysis type:	trib end point detection	
Date and time:	Fri, Jan 14 15:25:36 2011	
Virus		
Virus	Result	Controls
Herpes Simplex Virus 1	Negative	Passed
Herpes Simplex Virus 2	Negative	Passed
Varicella Zoster Virus	Negative	Passed
Epstein-Barr Virus	Negative	Passed

### ORDERING REFERENCES

- CLART® ENTERPEX Extraction  
48 tests: AT-0908-48
- CLART® ENTERPEX Amplification  
48 tests: AT-1008-48-MT
- CLART® ENTERPEX Visualization  
48 tests: CS-1108-48

> CLART® ENTERPEX: Essential features of CLART® ENTERPEX are protected by Patent Family of International PCT Patent application WO2009122201.

> CLART® ENTERPEX kit, accomplish with all the normatives described in EU 98/79/CE for IVD.


> CMV detection has been assented by the NB 0318.



# CLART® EnteroBac

## Detection of Bacteria causing Infectious Diarrhoea

### BACTERIA DETECTED



### MAIN ADVANTAGES OF MOLECULAR METHODS IN GASTROENTERIC DISEASES

- Direct detection of a wide spectrum of bacteria causing infectious diarrhoea.
- Detection of single or multiple infections of the most prevalent bacteria causing infectious diarrhoea in one test.
- Early detection of common pathogens causing outbreaks, such as *Salmonella*.
- No stool culture required.
- Classical methods such as culturing or biochemical and phenotypic tests may not detect coinfections.

### FEATURES

- The kit is validated for automatic DNA extraction directly from stool samples.
- High sensitivity and specificity.


- Three quality controls included per sample:
  - **Genomic DNA control:** validates the extraction performance.
  - **Amplification control:** avoids false negative results.
  - **Biotin markers:** check the proper performance of the visualization reagents provided with the kit.
- Each bacteria is detected in quadruplicate avoiding unspecific bindings.
- From sample to result in less than a working day.

### Data Management

- Automatic reading and interpretation of results within 5 hours after sampling.
- User-friendly report format (html, bmp).
- Samples are analyzed individually and three complementary reports are generated.
- Each report generated can be stored in the reader, exported or printed.

## REPORTING RESULTS

> Report and image obtained by CAR® reader.



B4		AT code: 60519	Rev. 10/13
Result view <b>CLART® EnteroBac</b>			
Sample reference:	3		
Array ID:	00000050516 (C1)		
Analysis type:	trib end point detection		
Date and time:	Fri, Jan 14 15:25:36 2011		
Bacteria	Result	Controls	
Aeromonas spp.	Negative	Passed	Passed
Shigella spp.	Positive	Passed	Passed
Salmonella spp.	Negative	Passed	Passed
Escherichia coli	Negative	Passed	Passed
Yersinia enterocolitica	Negative	Passed	Passed
Yersinia spp.	Negative	Passed	Passed
Campylobacter spp.	Negative	Passed	Passed
Clostridium difficile	Negative	Passed	Passed

## ORDERING REFERENCES

CLART® EnteroBac Amplification

48 tests: CS-0611-48

CLART® EnteroBac Visualization

48 tests: CS-0711-48


> CLART® EnteroBac kit, accomplish with all the normative described in EU 98/79/EC for IVD.



# CLART® SeptiBac

## Sepsis Detection

### BACTERIA AND FUNGI DETECTED



### MAIN ADVANTAGES OF MOLECULAR METHODS IN SEPSIS DIAGNOSTICS

Molecular methods reduce the limitations and drawbacks of conventional blood cultures:

- Blood/Sample volume.
- Turnaround time to definitive identification.
- Higher sensitivity to slow-growing and fastidious organisms.
- Co-infections detected
- Methicillin resistance marker included (*mecA* gene).

### FEATURES

- The kit is validated for automatic DNA extraction from positive blood cultures.
- High sensitivity and specificity.


- Three quality controls included per sample:
  - **Genomic DNA control:** validates the extraction performance.
  - **Amplification control:** avoids false negative results.
  - **Biotin markers:** check the proper performance of the visualization reagents provided with the kit.
- Targets are detected in triplicate avoiding unspecific bindings.
- Results obtained within a working day (4 hrs).
- Reduces turnaround time, allowing therapy adjustments.

### Data Management

- Automatic reading and interpretation of results.
- User-friendly report format (html, bmp).
- Samples are analyzed individually and three complementary reports are generated.
- Each report generated can be stored in the reader, exported or printed.

### REPORTING RESULTS

> Report and image obtained by CAR® reader.



A2		A7 min
Result view		0099
CLART® SeptiBac		Rev. 06/13
Sample reference:	3	
Array ID:	0000005016 (C1)	
Analysis type:	trib end point detection	
Date and time:	Fri, Jan 14 15:25:36 2011	
Bacteria	Result	Controls
<i>C. albicans</i>	Positive	Passed
<i>C. glabrata</i>	Negative	Passed
<i>C. Krusei</i>	Negative	Passed
<i>C. perfringens</i>	Negative	Passed

### ORDERING REFERENCES

CLART® SeptiBac Amplification

48 tests: CS-0311-48

CLART® SeptiBac Visualization

48 tests: CS-0411-48

> CLART® SeptiBac kit accomplish with all the normative described in EU 98/79/EC for IVD.

# BIBLIOGRAPHY

## CLART® CMA KRAS·BRAF·PI3K & NRAS·IKRAS

1. "Biomarkers Predicting Clinical Outcome of Epidermal Growth Factor Receptor-Targeted Therapy in Metastatic Colorectal Cancer". Review J Natl Cancer Inst 2009;101:1308-1324
2. "Effects of KRAS, BRAF, NRAS, and PIK3CA mutations on the efficacy of Cetuximab plus chemotherapy in chemotherapy-refractory metastatic colorectal cancer: a retrospective consortium analysis". Lancet Oncol 2010; 11: 753-62.
3. "KRAS, BRAF, PIK3CA, and PTEN mutations: implications for targeted therapies in metastatic colorectal cancer". Lancet Oncol 2011; 12: 594-603.
4. "Molecular Mechanisms of Resistance to Cetuximab and Panitumumab in Colorectal Cancer". J Clin Oncol 2010 28:1254-1261.
5. "NRAS and KRAS testing by a new diagnostic method to detect point mutations in colorectal cancer specimens: CLART NRAS-iKRAS. 2014 ASCO Annual Meeting Proceedings".

## CLART® HPV2

1. "Prevalence of HPV Infection among females in the United States". JAMA. 2007; 297(8):812-81.
2. "High frequency of multiple HPV types in cervical specimens from Danish women". APMIS 2009, 117: 108-114.
3. "External quality assessment for molecular detection of human papillomaviruses". Journal of Clinical Virology 48 (2010) 251-254.
4. "Human Papillomavirus 2 Assay Compared With the Hybrid Capture 2 Test". Journal of Medical Virology 2011 83:272-276 (2011).
5. "Identification of Multiple HPV Types on Spermatozoa from Human Sperm Donors". PLOS ONE. March 2011, Volume 6, Issue 3, e18095.

6. "Human papillomavirus genotype distribution among French women with and without cervical abnormalities". Intl J of Gynecol & Obstetrics 2011, Vol 114, Issue 2, Pag116-119.
7. "Prevalence of Human Papillomavirus Infection in Women in Portugal. The CLEOPATRE Portugal Study". Int J Gynecol Cancer 2011;21: 1150Y1158.
8. "Detection and genotype distribution of human papillomavirus (HPV) DNA in Danish colorectal carcinoma patients". Poster presented at the 28th IPVC, Puerto Rico, 2012.
9. "Human Papillomavirus Type Distribution in Cervical Intraepithelial Neoplasia Grade 2/3 and Cervical Cancer in Portugal. A CLEOPATRE II Study". Int J Gynecol Cancer 2013;23: 500Y506.

10. "Patterns of cervical coinfection with multiple human papilloma virus types in a screening population in Denmark". Vaccine, Volume 31, Issue 12, 15 March 2013, Pages 1604-1609.

## CLART® STIs A&B

1. "Gardnerella, Trichomonas vaginalis, Candida, Chlamydia trachomatis, Mycoplasma hominis and Ureaplasma urealyticum in the genital discharge of symptomatic fertile and asymptomatic infertile women". New Microbiologica, 33, 69-76, 2010.
2. "Global strategy for the prevention and control of sexually transmitted infections: 2006 -2015 : breaking the chain of transmission". WHO
3. "Sexually Transmitted Diseases in the United States, 2008 National Surveillance Data for Chlamydia, Gonorrhea, and Syphilis" CDC.
4. "Persistent increase in the incidence of acute male urethritis diagnosed in general practices in France". British Journal of General Practice 2006; 56: 110-114.
5. "Mycoplasma genitalium presence, resistance and epidemiology in Greenland". Int J Circumpolar Health 2012, 71:18203

## CLART® ENTHERPEX

1. "L-HHV7: co-facteur des infections méningées à entérovirus?". Poster presented at RICAI, December 2009.
2. "Evaluation of a new commercial PCR-DNA microarray for rapid and simultaneous detection of 9 viruses responsible for central nervous system infections". Poster presented at RICAI, December 2009.
3. "Combining Multiplex Reverse Transcription-PCR and a Diagnostic Microarray to Detect and Differentiate Enterovirus 71 and Coxsackievirus A16". Journal of Clinical Microbiology, 44(6), 2212-2219 (2006).
4. "DNA microarrays for virus detection in cases of central nervous system infection". Journal of Clinical Microbiology. 42(12), 5811-5818 (2004).
5. "Rapid Virological Diagnosis of Central Nervous System Infections by Use of a Multiplex Reverse Transcription-PCR DNA Microarray". Journal of Clinical Microbiology, Nov. 2011, p. 3874-3879 Vol. 49, No. 11.
6. "Human Herpesvirus-6A/B Binds to Spermatozoa Acrosome and Is the Most Prevalent Herpesvirus in Semen from Sperm Donors". PLOS ONE. November 2012, Vol. 7, Issue 11, e48810.

## CLART® PneumoVir

1. "Tract Viral infections and Coinfections in Patients with Influenza like Illnesses by use of RT-PCR DNA Microarray Systems". J. Clin. Microbiol. doi:10.1128/JCM.00733-10.
2. "Characterization of viruses causing Human Respiratory Infections via Genomic Identification for in vitro diagnosis. CLINICAL ARRAYS/CLART PneumoVir". Poster Session at ECCMID 2008, Barcelona.
3. "Myocarditis Caused by Human Parainfluenza Virus in an Immunocompetent Child Initially Associated with 2009 Influenza A (H1N1) Virus". Journal of Clinical Microbiology, May 2011, p. 2072-2073

4. "Respiratory viruses in Children Admitted to Hospital Intensive Care Units: Evaluating the CLART® PneumoVir DNA Array". Journal of Medical Virology 83: 150-155 (2011).
5. "Evaluation of viral co-infections in hospitalized and non-hospitalized children with respiratory infections using microarrays". Clin Microbiol Infect, 10.1111/1469-0691.12015.
6. "Co-infección viral respiratoria en niños hospitalizados por infección respiratoria aguda y su impacto en la gravedad clínica". Rev Chil Infect 2012; 29 (2): 169-174.
7. "Broad Respiratory Virus Detection in Infants Hospitalized for Bronchiolitis by Use of a Multiplex RT-PCR DNA Microarray System". Journal of Medical Virology 84:979-985 (2012)

## CLART® EnteroBac

1. "Diagnóstico microbiológico de infecciones gastrointestinales". 2nd Ed. (24). Procedimientos en Microbiología Clínica. Sociedad Española de Microbiología Clínica, SEIMC (2007).
2. "Diseño y optimización de un sistema rápido por multiplex-PCR e hibridación en microarray para la detección e identificación de patógenos bacterianos en pacientes con diarrea. Enfermedades Infecciosas y Microbiología Clínica. Vol 29, 33-34 (2011). 15th Congreso de la Sociedad Española de Microbiología Clínica (SEIMC), Junio 2011

## CLART® SeptiBac

1. "Hemocultivos". 2nd Ed. (3a). Eds. Cercenado E. and R. Canton. Procedimientos en Microbiología Clínica. Sociedad Española de Microbiología Clínica, SEIMC (2003).
2. "Diseño y optimización de un sistema de detección molecular por microarrays para la rápida identificación de bacterias gram positivas y hongos en hemocultivos positivos". Enfermedades Infecciosas y Microbiología Clínica. Vol 29, 33-34 (2011). 15th Congreso de la Sociedad Española de Microbiología Clínica (SEIMC), Junio 2011.