

Fully automatic system, completely guided by hybriSoft software











- √ Fully automated: AMPLIFICATION, DENATURATION and HYBRIDIZATION
- √ UV light for DNA decontamination
- √ Extremely robust and compact
- √ Faster, sensitive and more accurate in diagnosis.
- √ Thermostatic reaction chamber
- √ Vacuum pump control included
- √ Processed samples:
- Amplification from 1 to 24 samples per run
- Hybridization from 1 to 12 samples per run



√ HybriSoft: Specific software developed for hybriSpot platforms





HPV POSITIVE GENOTYPES DETECTED High Risk: 45

- Samples administration
- Automatic capture of images
- Automatic analysis of final arrays
- **Printed report**
- **Multitest**

Reference: VIT-HS12A

Developed and manufactured by



VITRO, S.A. Luis Fuentes Bejarano, 60 Edificio Nudo Norte (Local 3) 41020 Sevilla Telf: 902 366 974 www.vitro.bio

Distributed by